



EN Titanium Forceps

DE Titanium-Zange

ES Pinzas de titanio

Page 97 to 98





ATRAUMA

TITANIUM



Closing Pressure 220 +/- 10 gr
Length 30 mm
Jaw Length 9 mm
Weight 1 gr

DF-Ti-1038
Vena
Closing pressure
110 +/- 10 gr

Arteria



DF-Ti-1039
Vena
Closing pressure
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr

Length 35 mm

Jaw Length 13 mm

Weight 1 gr



DF-Ti-1040
Vena
Closing pressure
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr

Length 40 mm

Jaw Length 14 mm

Weight 1 gr



DF-Ti-1041
Vena
Closing pressure
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr

Length 50 mm

Jaw Length 23 mm

Weight 1.5 gr



Closing Pressure 220 +/- 10 gr
Length 30 mm
Jaw Length 9 mm
Weight 1 gr

DF-Ti-1042
Vena
Closing pressure
110 +/- 10 gr

Arteria



DF-Ti-1043
Vena
Closing pressure
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr

Length 35 mm

Jaw Length 13 mm

Weight 1 gr



DF-Ti-1044
Vena
Closing pressure
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr

Length 40 mm

Jaw Length 14 mm

Weight 1 gr



DF-Ti-1045
Vena
Closing pressure
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr

Length 50 mm

Jaw Length 23 mm

Weight 1.5 gr



Closing Pressure 220 +/- 10 gr
Length 30 mm
Jaw Length 9 mm
Weight 1 gr

DF-Ti-1046
Vena
Closing pressure
110 +/- 10 gr

Arteria



DF-Ti-1047
Vena
Closing pressure
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr

Length 35 mm

Jaw Length 13 mm

Weight 1 gr



DF-Ti-1048
Vena
Closing pressure
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr

Length 40 mm

Jaw Length 14 mm

Weight 1 gr



DF-Ti-1049
Vena
Closing pressure
110 +/- 10 gr

Arteria

Closing Pressure 220 +/- 10 gr

Length 50 mm


Jaw Length 23 mm

Weight 1.5 gr

ATRAUMA



DEBAKEY




DF-Ti-1001
16 cm/6 1/4"
1.5 mm




DF-Ti-1002
20,5 cm/8 1/8"
1.5 mm



DF-Ti-1003
24,5 cm/9 5/8"
1.5 mm



DF-Ti-1004
16 cm/6 1/4"
1.8 mm



DF-Ti-1005
20,5 cm/8 1/8"
1.8 mm



DF-Ti-1006
24,5 cm/9 5/8"
1.8 mm

